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USER'S GUIDE

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When using any electrical products, especially if children are present, certain basic safety precautions should always be followed, including the following:

DANGER To avoid the risk of electric shocks:

- Do not use this device in the bath or in the shower.
- Do not handle the USB charger with wet hands.
- Do not place the USB charger or drop into water or any other liquid.
- Check the condition of the USB charger power cord before first use, then regularly during the device life.

WARNING To prevent against the risk of fire, burn, or injury to persons:

- Do not connect this device to a voltage system other than the one specified on the device.
- Do not use this device if you notice that a cord or a plug is damaged.
- Do not direct the water under the tongue, inside the ear, nostrils or other delicate areas. This device is capable of producing pressure that could cause serious injury in these areas.
- Use the device as indicated in the instructions for use or follow the recommendations of your dentist, periodontist or pharmacist.
- Use only tips and accessories recommended by NEOPULSE®.
- Do not drop the device
- Do not insert objects into the device holes.
- Keep the USB charger away from heated surfaces.
- Do not use outdoors or operate the device, where aerosol products are being used or where oxygen is being administered.
- Do not use concentrated insoluble essential oil in this device. The use of these products will degrade product performance and longevity. The use of these products may cancel the warranty.
- Remove any oral jewelry items before use.
- Do not use the device if you have an open wound on the tongue, in the mouth or in the nose.
- If your doctor has advised you to take antibiotics before dental care, consult your dentist, periodontist or pharmacist before using this device or any other oral hygiene aid.
- Give specific instructions to children (from 8 years) and people with special needs for proper use of this device and supervise them.
- This device does not contain any userreplaceable components, and does not require maintenance.
- Fill the tank only with lukewarm water or another solution recommended by a healthcare professional. Be careful that it does not contain dyes that could damage the system.
- If the product is opened or taken apart for any reason, the warranty will be invalidated.

2 YEARS WARRANTY

NEOPULSE warrants that this product is free from material and manufacturing defects for a period of two years from the date of purchase.

Please contact NEOPULSE: conso@neopulse.fr

PRODUCT DESCRIPTION

1. Tip
2. Power on/off button
3. Pressure control selection button
4. LED indicators, 6 pressure levels
5. Battery charge indicator
6. Charger plug
7. Removable tank
8. Tank filling lid

Accessoires :

- USB charging cable
- Classic tip (x2): Perfect for daily use (to replace every 6 months)
- Orthodontic tip (x2): Perfect to clean orthodontic appliance (braces) and for general use (to replace every 3 months)

Pressure settings:

You can adjust the pressure on the dental water jet by successive presses on the pressure control button (3). The LED indicators (4) shows the pressure level. 6 pressure levels: 2 levels for Normal, 2 levels for Soft and 2 levels for Pulse

Normal pressure: ideal for those who are used to dental water jet and who are looking for power for better efficiency:

- Normal 1 (7.6 bars): more intense blue light
- Normal 2 (6.1 bars): less intense blue light

Soft pressure: ideal for novices to get used to the pressure of the water, for children (+8 years) and people with sensitive teeth:

- Soft 1 (4.75 bars): more intense blue light
- Soft 2 (1.7 bars): less intense blue light

Pulse pressure: ideal for gum massage:

- Pulse 1 (6.9 bars): more intense blue light
- Pulse 2 (4.13 bars): less intense blue light

Recommendation: Start with the lowest setting and increase the pressure according to your own preferences.